



Labels



Freezerworks
Sample Management Software

Wrap-Around Labels



Environmental Conditions:

Resistance to:

- Xylene
- Alcohol
- Liquid Nitrogen
- Formalyn


Temperature Range:

-196°C to 149°C

Material:

- Vinyl label flood coated on half for thermal transfer printing.
- Frosty clear material wrapped around for a self laminating hold that protects the white printed area after application. Great for long term storage.


Sizes:



1-2 ml vial size Portrait Wraparound (self laminating)


LBL-17WR: 2.625" x 1.00" (66.675mm x 25.4mm) Printable area is 1.00" x 1.00" (25.4mm x 25.4mm). Fits most cryovials.

Packaged for printers with 3" label core roll (e.g., Cab EOS1, Zebra ZM400 printers)



1-2 ml vial size Portrait Wraparound (self laminating)


LBL-02WR: 2.75" x 1.00" (69.85mm x 25.4mm). Printable area is 1.3" x 1.00" (33.02mm x 25.4mm). Fits most cryovials. Use ribbon RIB-04 or RIB-05 with desktop printers (one ribbon per approx 1.5 roll of 600).



1-2 ml vial size Landscape


LBL-13WR: 2.75" x 1.0625" (70mm x 27mm). Printable area is 1.5" x 1.06" (38mm x 27mm).

Fits most cryovials. Use RIB-05 Ribbon for desktop printers. (one ribbon per approx. 3 rolls of 800).



.5 ml vial size

LBL-09WR: .625" x 2.25" (15.88mm x 57.15mm). Printable area 0.625" x 0.75" (15.88mm x 19mm). Portrait label that fits smaller, conical vials. For desktop printers with 1" (25.4mm) label core diameter. Use ribbon RIB-04 or RIB-05 (one ribbon per approx. 1.75 roll of 600).



LBL-19WR: .625" x 2.375" (15.88mm x 60.325mm) .5ml labels rolled landscape with 3" core for larger printers (e.g., cab EOS1, Zebra ZM400). Printable area .625" x 1.00" (15.88mm x 25.4mm). 2,000 per roll.

Other Label Formats:



Combo labels: Use for labelling the tops and sides of vials.

LBL-06COM, or LBL-07FBCOM for Freezerbondz material. Available in Freezerbondz™ material to adhere to frozen vials.

LBL-07FBS – Freezerbondz combo label on 1" (25.4mm) label core for smaller desktop printers

Good for -85C freezer storage

Dimensions: vial side: 26mm L x .375mm H + 9.530mm vial cap label diameter; 0.9" L x 0.9" H + 0.375" cap label diameter

LBL-18COM1 (1" or 25.4mm label core); **LBL-18COM3** (3" or 76.2mm label core) – **NitroTAG™**, performs best in **LN2 (-196C)** Dimensions: vial side: 23.9mm L x 12.7mm H + 11mm vial cap diameter; 0.94" L x 0.5" H + 0.437"



Freezerbondz™ Labels

- 1" x 1" (25.4mm x 25.4mm) Freezerbondz™: Labels adhere to frozen vials, lowering the risk of thaw when labeling samples.
- **LBL-05FBL** (3" or 76.2mm label core diameter for larger printers).
- **LBL-05FBS** (1" or 25.4mm label core for smaller desktop printers)
- **LBL-16:** 1.625" x .6" (41.275mm L x 15.24mm H) Freezerbondz™ Label: Adhere to frozen vials.



Labels for Slides:

- **LBL-15:** An .875"x.95" (22.225mm x 24.13mm) label with a clear plastic overcoat "flap" for slides. Withstands the chemical process and harsh environment of pathology staining.
- **LBL-08SL:** 0.9" x 0.9" Slide label Part #: LBL-08SL.

Paper Labels:

- **LBL-14** 2.25" x 1.25" (57.15mm x 31.75mm). A paper label for non-cryogenic uses. Roll of 2100.



Ribbons:

- **RIB-05:** 4.375" x 242' (111mm x 74M) wax/resin roll for desktop printers (e.g., Zebra GX420t, Brady BBP-11, TSC TTP 343C Printers).
- **RIB-04:** 2.52" X 242' (64mm x 74M) wax/resin roll for desktop printers (e.g., Zebra GX420t, Brady BBP-11, TSC TTP 343C printers). For printing all labels except landscape wraparound (LBL-13WR) labels.
- **RIB-06:** 3.15" x 1476' (80mm x 450M) wax/resin roll for larger printers (e.g., cab EOS1, Zebra ZM400) for landscape wraparound.
- **RIB-4311:** 1" x 984' (25.4mm x 300M) wax/resin roll for larger printers (e.g., ZM400, cab EOS1) for Portrait wraparound, 1'x1' Freezerbondz, Combo and Combo Freezerbondz.
- **RIB-6400:** 2.360" x 984' (60mm x 300M) resin roll for larger printers (e.g., cab EOS1, ZM400) for 1'x1' Freezerbondz, Combo and Combo Freezerbondz.